

Digital Neonatal Scales



Healthweigh™

H611-00-1, H610-00-1; H611-002, H610-002; H611-004, H610-004

Digital Neonatal Scale

Model	Capacity	Graduation	Product Weight	Product Dimensions
Class III	12kg 26lb	±5g ±0.2oz	9.5kg 21lb	340 x 560 x 200mm 13" x 22" x 7.75"
H611-00-1				
H611-00-2				
H611-00-4				



H611-00-1

Healthweigh™ Neonatal scale is designed with advanced technology for precision weighing of newborn and growing babies. The scale is equipped with movement compensation technology, essential to eliminate involuntary movement caused by baby. Designed for multifunctional tasking, the scale tray has a built-in measuring tape to measure and weigh baby at the same time. Lift-off and removable tray for easy cleaning, the tray comes in three options:

H611-00-1 opaque white plastic

H611-00-2 two side crystal clear plastic

H611-00-4 four side crystal clear plastic

Optional for user convenience is a 12 drawer durable plastic trolley on roller wheels, custom designed for the Healthweigh™ Neonatal scale.



H611-00-1; H610-00-1

Features

Movement compensation technology

Reweigh;

Zero

Hold

kg/lb

Healthweigh™ Sterling 6 digit 1" LCD

Updatable firmware

Power

AC230V/120V Input DC9V Output

Auto Power Off

Rechargeable battery / Low battery indication

Communication

Optional: RS232 port

NEW USB port; Bluetooth for wireless



H611-00-2; H610-00-2

Non-class III

Model	Capacity	Graduation
H610-00-1	0-2kg	±1g
H610-00-2	0-5lb	±0.05oz
H610-00-4	2-6kg	±2g
	5-12lb	±0.1oz
	6-15kg	±5g
	12-30lb	±0.2oz



H611-00-4; H610-00-2