

Digital Veterinary Scale



Healthweigh™ H741-00-2, H740-00-2; H711-00-2, H710-00-2 Digital Veterinary Scale

Model	Capacity	Graduation	Product Weight	Product Dimensions
Class III	220kg	100g	20kg	1066 x 560 x 50mm
H741-00-2	485lb	0.2lb	44lb	42" x 22" x 2"
H711-00-2				

The Healthweigh™ digital Veterinary scale is a versatile scale for weighing medium to large size animals. The scales are equipped with movement compensation technology, essential to eliminate involuntary movement caused during the weigh process.

The metal platform is protected with epoxy coating and guaranteed non-slip by an ultra fine, durable rubber mat. For convenience and hygiene, the rubber mat is removable, and extremely easy to clean.

The digital Veterinary Scales comes with Healthweigh™ Platinum Indicator, with tare function to predetermine weight of collar, basket or other item accompanying pet.

Features

Movement compensation technology
Reweigh;
On-Off
Zero
Hold
Tare
kg/lb
Updatable firmware

Optional

Healthweigh™ Sterling indicator 6 digit 1" LCD

Power

AC230V/120V Input DC9V Output
Auto Power Off
Rechargeable battery / Low battery indication

Communication

Optional: RS232 port
NEW USB port; Bluetooth for wireless

Non-class III

Model	Capacity	Graduation
H740-00-2	310kg	100g
H710-00-2	700lb	0.2lb



H741-00-2, H740-00-2