



Healthweigh™ Digital Indicator

Healthweigh™ Platinum Digital Indicator

The Healthweigh™ Platinum Digital Indicator is designed for Shekel Healthweigh™ health and fitness scales. The Platinum Digital Indicator is the standard indicator for most of the Healthweigh™ digital scale collection. The 6 digit 1" LCD display has large easy to read digits.

Four large indicator keys for:

- On-Off/Zero
- BMI (body mass index) calculation / Tare
- Reweigh / Scroll Down
- kg/lb / Scroll Up

Features

Low battery indication / automatic power-off
Hold function
AC Adaptor plug-in storage is an innovative feature designed for enhanced user convenience.

Power

AC230/120V Input; DC9V Output
Rechargeable battery

Communication

RS232 port enables data transmission for full integration into digital medical systems
NEW: Bluetooth for wireless application
USB port for data transmission



Platinum Indicator



AC adaptor plug-in storage

Healthweigh™ Sterling Digital Indicator

The Healthweigh™ Sterling Digital Indicator is designed for Shekel Healthweigh™ health and fitness scales. It is easily integrated into any scale or used as stand-alone on bracket for table or wall mount. The slim line design 6 digit 1" LCD is easy to read.

The four large keys are simple to use:

- On-Off
- Zero
- Reweigh
- kg/lb selection

Features

Low battery indication / automatic power-off
Hold function

Communication

Optional: RS232 port enables data transmission for full integration into digital medical systems
NEW: Bluetooth for wireless application
USB port for data transmission

Power

AC230/120V Input; DC9V Output
Rechargeable battery



Sterling Indicator